



AIRWAY

---

## Portrait Style No. 1650

---

**P**ortrait incorporates a totally new concept of weight distribution in a symmetrically shaped breast form. Its unique design properly distributes the weight from the chest wall flowing to the apex enabling the form to react naturally with body movements.

1650 BK now available to meet the needs of the darker skin tones.

---



---

## TriLite Style No. 1654

---

**T**he TriLite is a truly new concept in a lightweight prosthesis. Its unique posterior surface design makes it especially suited for wear while swimming, as it permits unimpaired drainage upon leaving the water and thus allows for rapid drying of skin. It provides perfectly natural contours, especially in a wet or dry bathing suit, and in lightweight lingerie. Weighs less than a comparable silicone breast form.

---



A naturally soft, resilient but stable shape with a perfect fit to the chest. Three silicone ribs reinforce the shell, creating a novel air chamber system that ensures rapid drying and optimum air circulation. Welcome weight relief, especially in cases of lymphadema or heavy breasts.

Size range accommodates 32A through 46C. 24 month limited warranty.

---

---

**TriShell** Style No. 1656

---

**T**riShell is comfortable, lightweight and conforms to the body and the bra for a totally natural feeling. It is ideally used as a lightweight form, for augmentations, reconstruction adjustments and for women suffering from congenital anomalies where one breast has not developed completely. Its three piece construction allows it to be used as a total breast replacement, as a partial, or as an enhancer by using only the soft silicone shell.

---



Soft outer silicone shell looks and moves like the natural breast. Inner cotton filler fashions itself with the chest wall creating a soft backing against the body. Double thick cotton back cover, molded seamless front.

Size range accommodates 32A through 46 D.

No. 1681: cover for TriShell, sizes 3-12

No. 1658: filler for TriShell, sizes 3-12

No. 1656: silicone shell only

---

