

# JAZZY ELECTRONICS

## JOYSTICK CONTROLLERS & OPTIONS



Remote Joystick Controller



Pilot Joystick Controller



Dynamic Europa Controller



Dynamic Dolphin Controller

Pride offers a variety of alternative joysticks knobs for your Jazzy PG8 controller system. In addition to the endless variety of chair/seating configurations necessary to fit patients' needs, special provisions may be required to operate joystick controls. Depending on the individual's special needs, the assortment of alternative knobs offered should prove to be a resource worth using. Use of the correct knob design for your customer will reduce fatigue and increase comfort to everyday driving. All knobs are designed to be installed with minimal effort and are 100% compatible with all Jazzy PG8 controllers.



Extended I-Handle



T-Handle



Mushroom Handle



Medium Ball Handle



Standard Handle

### Swing-Away Joystick

The Swing-Away mount allows for the remote joystick controller to be locked in either a forward mount or a retracted position allowing closer access to tables and desks. This unique design allows for the remote joystick to remain oriented in a normal drive mode configuration. The swing away mount is compatible with all seating options and is fully retrofittable to all previous remote mounts.

