

JAZZY SEATING OPTIONS

MEDIUM-BACK & HIGH-BACK SEATING

Medium-Back Seat

The Medium-Back style seat is designed to be used for patients that are not at risk for skin breakdown and have little or no special positioning needs. This seating option makes the transition from a scooter to a power chair much easier for these patients.



*Jazzy 1100 with
Medium-Back Seat
in Champagne.*



*Jazzy 1100 with Gray Deluxe
High-Back Seat in Red.*

Deluxe High-Back Seat

The Deluxe High-Back Seat offers greater comfort than the standard Medium-Back seat due to the higher back with lumbar support and deeper contoured seat cushion. The High-Back Seat is designed for patients not at risk for skin break down or patients that don't have highly involved positioning issues. This seating option makes the transition from a scooter to a power chair much easier for these patients.

New Solid Seat Pan for Medium-Back and High-Back Seats

The Solid Seat Pan offers the Medium and High-Back seats the ability to accommodate clients with positioning needs and are at risk of skin break down. The Solid Seat Pan is available in various sizes and is covered in black vinyl. Black velcro strips are provided for easy mounting of seating systems. Back cushions are available in all standard colors for Medium-Back and High-Back seats.



*Jazzy 1100 with Beige Deluxe
High-Back Solid Seat Pan and
ROHOTM cushion in
Caribbean Teal*

JAZZY SEATING OPTIONS

POWER SEAT & RECLINING SEAT



*Jazzy 1100 with Beige
Deluxe High-Back Seat in
Forest Green.*

Power Seat

The Jazzy Power Seat gives the user the convenience of improved accessibility to his or her daily environment. The Power Seat option is adjustable from 19" to 25" and is available with all seating configurations. This provides a total of 6" of seat lift over the standard height, allowing easy access to objects not reachable with a standard fixed seat. The Jazzy Power Seat option can be controlled by either a toggle switch or an optional accessory remote joystick. This option also allows for a seat swivel of 90 degrees providing for easier transfers.

Reclining Seat

The Jazzy Reclining Seat is adjustable from a standard 90 degrees to a recline position of 180 degrees. The Recline mechanism is designed to be infinitely adjustable, providing customized comfort. The Reclining Seat is also available with a solid seat pan to use with pressure relieving cushions.



*Jazzy 1120 with Gray
Reclining Seat in
Viper Blue.*

*Jazzy Model 1105,
shown with
Specialty Seat,
Pilot Controller in
Arctic Blue.*



Specialty Seating

The Specialty Seat is an economy seat frame that offers adjustable seat width of 16", 18" and 20". The Specialty Seat is designed with 7/8" tubing and will be able to accept all types of seating systems. The Specialty seat comes standard with adjustable height and angle flip back arm rests.

JAZZY SEATING OPTIONS

VERSA SEATING SYSTEMS

Versa Seat

The most innovative and adjustable seating system available in the industry today.

Pride's Versa seat is the most innovative and adjustable seating system on any Power Chair available today. This Rehab seating system provides you, the rehab professional, with the ability to make a full range of height, width, depth and angle adjustments "on site" according to the exact requirements of the end user. In addition, the Versa Seat can be periodically readjusted to change with the needs of the end user. The Versa Seat interfaces with all types of seating systems making it the most user friendly system on the market today. The innovative Versa Seat is available exclusively on the #1 rated Jazzy Power Chair.



Jazzy 1100 with Versa Seat and Swing-Away Footrests



Design Features Include:

- Removable Solid Seat Base Pan that Provides Stable Surface for Seating Interfaces
- Accepts Pride's Swing Away Footrests and Elevating Leg Rests
- Seat Base Adjustable in Seat Width (13" to 21" in 1" increments) and Seat Depth (10" to 24" in 1" increments.)
- Angle Adjustable Back with 24 Degrees of total Adjustment
- Height and Angle Adjustable Seat Pan with Adjustable Center of Gravity with 5" of height adjustability and 20 degrees of angle adjustability
- Interfaces with all Industry Seating Systems
- Black Powder Coat Finish Offers Attractive Styling and Durability
- Retro-Fittable to the Complete Jazzy Power Chair Line
- Works with Most Jazzy Power Chair Options
- Angle and Height Adjustable Flip-up Arms
- Detachable, Height Adjustable Arms
- Durable 1" Square Tube Frame Construction Provides Strength and Adjustability
- Rear Canes Constructed of 1" Tubing to Interface with all Seating Systems

VERSA SEATING SYSTEM ARMRESTS



Flip back height adjustable. Standard armrest available in desk length and full length.

JAZZY SEATING OPTIONS

VERSA SEATING SYSTEMS



New Versa Tilt Seating Systems

Introducing Pride's new Versa Tilt CG seating system for Jazzy Power Chairs. The Versa-Tilt CG system changes the center of gravity as the client tilts backwards. This allows the Jazzy to maintain the best turning radius in the industry. Since the Versa-Tilt CG is retrofittable to fit the Jazzy 1100, Jazzy 1120 and the Jazzy XL 1170, it will allow you to transition from any Jazzy seating option to a power tilt seating system while keeping the same power base. The added advantage of the Versa-Tilt system is that it utilizes Pride's Versa Seat technology which allows the seating system to be adjustable in the field without the need to purchase any additional parts. The solid base has an adjustable width of 16" – 22" and an adjustable depth of 14" – 22" in 1" increments.

Design Features Include:

- Weight Shifting C.G. Tilt (0 – 60 degrees or 5 – 60 degrees)
- Removable Solid Seat Base Pan that Provides Stable Surface for Seating Interfaces
- Accepts Pride's Swing Away Footrests and Elevating Leg Rests
- Seat Base Adjustable in Seat Width (16" to 22" in 1" increments) and Seat Depth (14" to 22" in 1" increments.)
- Angle Adjustable Back with 24 Degrees of total Adjustment
- Seat Height Adjustable up to 4" and Seat Angle Adjustable from 0 – 5 degrees fixed tilt.
- Interfaces with all Industry Seating Systems
- Black Powder Coat Finish Offers Attractive Styling and Durability
- Works with Most Jazzy Power Chair Options
- Angle and Height Adjustable Flip-up Arms or Detachable Arms
- Durable 1" Square Tube Frame Construction Provides Strength and Adjustability
- Rear Canes Constructed of 1" Tubing to Interface with all Seating Systems
- Low Current Switching and Drive Disable Safety Switch
- 5 Year Frame Warranty
- 18 Months Electronics and Actuator Warranty

VERSA SEATING SYSTEM ARMRESTS



Removeable, height adjustable.

Optional armrest available in desk length and full length.

**Custom
Painting
Available!**



Jazzy Model 1170 X with Versa Seat and Swing-Away Footrests (Custom Painted Shroud)