

# JAZZY 1105



*Jazzy 1105 shown in  
Arctic Blue with Gray  
High Back Seat and  
Swing Away Footrests*



The Jazzy model 1105 is a stylish, high performance foldable Power Chair which is easy to transport. The Jazzy 1105's mid-wheel drive design provides a smooth and comfortable ride, superb maneuverability and outstanding handling and control.

***Jazzy Power Chairs®***  
***How a Power Chair Should Feel™***

## SPECIFICATIONS

Wheels:	14" pneumatic, center mounted
Rear Casters:	8" solid
Anti-tip:	6" solid, front mounted
Maximum Speed:	4.5 MPH
Ground Clearance:	3.5 inches
Turning Radius:	19.5 inches
Overall Size:	Length: 39.5 inches Width: 23 inches
Seating Options:	Medium Back (standard), Medium Back with Solid Seat Pan Width: 17" Depth: 16" High Back, High Back with Solid Seat Pan Width: 17" Depth: 18" Manual Recline, Versa Seat
Drive Train:	Two motor mid-wheel drive
Battery:	Two 12 volt, 22NF batteries
Range:	Up to 25 miles
Battery Charger:	On-board (standard)
Electronics:	50 amp pilot controller
Weight Capacity:	250 lbs.
Transport Weight:	Heaviest piece 71 lbs.
Color:	Evergreen or Arctic Blue
Warranty:	Lifetime limited warranty on frame Two year electronics warranty 18 month on drive motors



Easy to transport.  
It can be carried  
or rolled.



## FEATURES

- Available in Evergreen or Arctic Blue
- Powerful two motor, mid-wheel drive design
- Patent pending mid-wheel drive design provides *outstanding* maneuverability
- Low center of gravity design provides a very smooth and stable ride
- 14" knobby mid drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4.5 m.p.h.
- Adjustable seat height and angle
- On-board battery charger
- Cam-Strut foldable frame
- 23 inch width makes it ideal for indoor and tight quarter maneuvering



Disassembles in four pieces

## OPTIONS

- Lap Belt
- Cupholder
- Dust Cover
- Rear Basket
- Oxygen Tank Holder
- Cane/Crutch Holder
- Walker Holder
- Flat Free Inserts
- Swing Away Leg Rests
- Elevating Leg Rests
- Heel Loops
- Adjustable Ht. Armrests

**Pride**  
Health Care, Inc

182 Susquehanna Ave. Exeter, PA 18643  
1-800-800-8586 1-800-800-1636 (fax)  
[www.PRIDEHEALTH.com](http://www.PRIDEHEALTH.com)

