

# OPTIONS & ACCESSORIES

# OPTIONS & ACCESSORIES

## Shroud Color Options

The Jazzy's aerodynamically contoured shroud is molded of the finest ABS plastic to provide high-impact strength. Coated in a durable, resilient Durathane DU-100 paint, the Pride Jazzy is available in a beautiful array of color



Candy Apple Red    Caribbean Teal    Forest Green    Gray    Purple Pizazz    Champagne    Viper Blue

## Footriggering Options

The Jazzy is available with a standard Foot Platform or optional Swing Away Foot Rest and Swing Away Elevating Leg Rest. These Foot Riggerings offer a high level of durability and adjustability. They are available on all Jazzy seating systems. All Footriggering options are finished in a black powder coat paint for a more durable finish and to compliment the Jazzy's elegant styling. They also feature a new camlock design for dependability and ease of use.



Elevating Leg Rest    Standard Foot Platform    Swing Away Footrest

## JAZZY SPECIFICATIONS

	<i>Jazzy 1103</i>	<i>Jazzy 1104</i>	<i>Jazzy 1115</i>	<i>Jazzy 1105</i>	<i>Jazzy 1100</i>	<i>Jazzy 1120</i>	<i>Jazzy 1170XL</i>
<b>Primary Use</b>	For a client that requires superb indoor maneuverability and limited outdoor use	For a client that wants elegant style and good maneuverability at an economical price	For a client that wants great indoor/outdoor and requires a transportable Power Chair	For a client that wants great indoor/outdoor and requires a transportable Power Chair	For a client with rehab needs that require excellent indoor and outdoor maneuverability	For the active rehab client providing the best balance of indoor and outdoor maneuverability	For a highly active rehab client requiring superb outdoor performance with increased weight capacity
<b>Wheels</b>	10" pneumatic	14" pneumatic	14" pneumatic	14" pneumatic	14" pneumatic	14" pneumatic	16" pneumatic
<b>Rear Casters</b>	8" solid	8" solid	8" solid	8" solid	8" solid	8" solid	9" pneumatic
<b>Anti-Tip Wheels</b>	6" solid	6" solid	6" solid	6" solid	6" solid	6" solid	8" solid
<b>Maximum Speed</b>	4 MPH	4.5 MPH	4.5 MPH	4.5 MPH	4.5 MPH	4.5 MPH or 6 MPH <sup>5</sup>	6 MPH
<b>Ground Clearance</b>	3 inches	3.5 inches	3.5 inches	3.5 inches	3.5 inches	3.5 inches	4.25 inches
<b>Turning Radius</b>	17 inches	19.5 inches	19.5 inches	19.5 inches	19.5 inches	19.5 inches	22.5 inches
<b>Overall Length</b>	34 inches	39.5 inches	39.5 inches	39.5 inches	39.5 inches	39.5 inches	45 inches
<b>Overall Width</b>	23 inches	25.4 inches	23 inches	23 inches	25.4 inches	25.4 inches	26.5 inches
<b>Medium Back Seat</b>	Width: 18" Depth: 16"	Width: 18" Depth: 16"	Width: 18" Depth: 16"	Width: 18" Depth: 16"	Width: 18" Depth: 16"	Width: 18" Depth: 16"	Width: 18" Depth: 16"
<b>Deluxe High Back Seat</b>	Width: 18" Depth: 16" <sup>2</sup>	Width: 18" Depth: 16" <sup>2</sup>	Width: 18" Depth: 16" <sup>2</sup>	Width: 18" Depth: 16" <sup>2</sup>	Width: 18" Depth: 16" <sup>2</sup>	Width: 18" Depth: 16" <sup>2</sup>	Width: 18" Depth: 16" <sup>2</sup>
<b>Versa Seat</b>	N/A	Width: 13"-21" Depth: 10"-21"	Width: 13"-21" Depth: 10"-21"	Width: 13"-21" Depth: 10"-21"	Width: 13"-21" Depth: 10"-24"	Width: 13"-21" Depth: 10"-24" <sup>3</sup>	Width: 13"-21" Depth: 10"-24" <sup>3</sup>
<b>Drive Train</b>	2-motor mid-wheel drive	2-motor mid-wheel drive	2-motor mid-wheel drive	2-motor mid-wheel drive	2-motor mid-wheel drive	2-motor mid-wheel drive	2-motor mid-wheel drive
<b>Battery</b>	Two 12-volt, U-1	Two 12-volt, group 24	Two 12-volt, 22NF	Two 12-volt, 22NF	Two 12-volt, group 24	Two 12-volt, group 24	Two 12-volt, group 24
<b>Per-charge Range</b>	Up to 25 miles	Up to 25 miles	Up to 25 miles	Up to 25 miles	Up to 25 miles	Up to 25 miles	Up to 25 miles
<b>Electronics</b>	50 amp Pilot controller	50 amp Pilot controller	50 amp Pilot controller	50 amp Pilot controller	Remote 70 amp controller or 50 amp Pilot controller	Remote 70 amp controller or 50 amp Pilot controller	70 amp Dynamic or P&G controller
<b>Weight Capacity</b>	250 lbs.	250 lbs.	250 lbs.	250 lbs.	300 lbs.	300 lbs.	400 lbs.
<b>Base Weight</b>	71 lbs.	81 lbs.	71 lbs.	71 lbs.	81 lbs.	81 lbs.	124.5 lbs.
<b>Seat Weight</b>	29 lbs.	29 lbs.	29 lbs.	27 lbs.	29 lbs.	29 lbs.	29 lbs.
<b>Battery Weight</b>	24 lbs. (each)	48 lbs. (each)	38 lbs. (each)	38 lbs. (each)	48 lbs. (each)	48 lbs. (each)	48 lbs. (each)
<b>Full Suspension</b>	No	No	No	No	Yes	Yes (Active-Trac) <sup>1</sup>	Yes (Active-Trac) <sup>1</sup>
<b>Curb Climbing Height</b>	1"	1"	1"	1"	2"	4"	6"
<b>Specialty Controls</b>	No	No	No	No	Yes w/ Remote Controller	Yes w/ Remote Controller	Yes
<b>Tilt in Space &amp; Recline</b>	No	No	No	No	Tilt & Recline	Tilt & Recline	Tilt & Recline

All models come with a lifetime limited warranty on frame, two year limited warranty on electronics and 18 month limited warranty on drive motor. Also standard is an on-board battery charger.

- <sup>1</sup> - The Active-Trac is a unique suspension system that allows the Jazzy to better traverse uneven terrain and ramp transitions.
- <sup>2</sup> - Standard size shown. Custom sizes available.
- <sup>3</sup> - Available up to 24 inches in width.
- <sup>4</sup> - With Heavy Duty Package
- <sup>5</sup> - 6 MPH top speed available with high speed motor package on 1120

## OPTIONS

- Available in seven stylish colors
- Deluxe High-Back Seat
- Power Seat
- Reclining Seat
- Versa Seat
- Specialty Seating
- Adjustable Height Arms
- Flat Free Tire Inserts
- Seat Belt
- Cupholder (L & XL)
- Rear Basket
- Basic Lighting Package
- Full Lighting Package
- Heavy Duty Dust Cover
- Swing Away Joystick
- Swing Away Footrest
- Elevating Leg Rest
- Heel Loops
- Oxygen Tank Holder
- Crutch & Cane Holder
- Walker Holder