

# SPECIFICATIONS



	CELEBRITY	LEGEND	XL SERIES	SUNDANCER	SIDEKICK
<b>Models</b>	SC-400, SC-401, SC-402 SC-440, SC-441, SC-442	SC-300, SC-340	SC-345, SC-445 SC-446, SC-447	SC-202, SC-203, SC-205	SC-200, SC-201
<b>Overall Size</b>	Length: 46" Width: 24"	Length: 46" Width: 24"	Length: 52" Width: 24"	Length: 41" Width: 22"	Length: 41" Width: 22"
<b>Range</b>	20-25 miles per battery charge	20-25 miles per battery charge	25-30 (with 40AH Battery)	20-25 miles per battery charge	20-25 miles per battery charge
<b>Speed</b>	5.25 mph (60% reverse)	5.25 mph (60% reverse)	5.7 mph (60% reverse)	4.5 mph (60% reverse)	4.5 mph (60% reverse)
<b>Turning Radius</b>	34"	34"	36"	31"	31"
<b>(w/o batteries) Weight</b> (31 AH Battery: 24 lbs.) (40 AH Battery: 31 lbs.)	Rear: 39.8 lbs Seat: 23 lbs. 3-Wheel Base/Tiller: 32.5 lbs. 4-Wheel Base/Tiller: 38.5 lbs.	Rear: 39.8 lbs. Seat: 23 lbs. 3-Wheel Base/Tiller: 32.5 lbs. 4-Wheel Base/Tiller: 38.5 lbs.	Rear: 49.5 lbs. Seat: 37 lbs. Base: 52.5 lbs.	Rear: 38 lbs. Seat: 23 lbs. Base: 27 lbs.	Rear: 38 lbs. Seat: 23 lbs. Base: 27 lbs.
<b>Ground Clearance</b>	4.0"	4.0"	5.5"	3.5"	3.5"
<b>Seat to Ground</b>	21.5" - 25.5"	21.5" - 25.5"	25" - 29"	18" - 20"	18" - 20"
<b>Front Tires</b>	(Pneumatic) 3" x 9"	(Pneumatic) 3" x 9"	(Pneumatic) 4" x 11.5"	(Pneumatic) 2" x 8"	(Pneumatic) 2" x 8"
<b>Rear Tires</b>	4" x 10"	4" x 10"	4" x 12"	2.5" x 9"	2.5" x 9"
<b>Battery</b>	Two 12-Volt Sealed	Two 12-Volt Sealed	Two 12-Volt Sealed	Two 12-Volt Sealed	Two 12-Volt Sealed
<b>Battery Charger</b>	On-board, 3.5 amp	On-board, 3.5 amp	On-board, 3.5 amp	On-board, 3.5 amp	On-board, 3.5 amp
<b>Rear Wheel Drive</b>	24-Volt DC Sealed Transaxle Motor	24-Volt DC Transaxle Motor	24-Volt DC Sealed Transaxle Motor	24-Volt DC Sealed Transaxle Motor	24-Volt DC Sealed Transaxle Motor
<b>Weight Capacity</b>	300 lbs. (400 lbs. Optional)	300 lbs. (400 lbs. Optional)	300 lbs. (400 lbs. Optional)	300 lbs. (400 lbs. Optional)	300 lbs. (400 lbs. Optional)
<b>Seat Colors</b>	Gray or Beige in Fabric or Vinyl Black Vinyl	Vinyl: Beige, Burgundy or Black Fabric: Fawn or Burgundy	Legend XL: Vinyl or Fabric in Burgundy or Fawn Celebrity XL: Vinyl or Fabric in Charcoal Gray Limited Edition XL: Vinyl in Beige, Fabric in Gold	Gray Fabric, Gray Vinyl, Black Vinyl	Vinyl: Brown, Burgundy, Black Fabric: Burgundy, Fawn
<b>Body Colors</b>	Viper Blue, Forest Green, Candy Apple Red	Champagne	Candy Apple Red, Forest Green, Viper Blue, Champagne (Legend)	Viper Blue, Teal, Candy Apple Red	Gray or Champagne
<b>Additional Features</b>	User Adjustable Sliding Seat System		User Adjustable Sliding Seat System		
<b>Disassembly</b>	5 easy to transport pieces	5 easy to transport pieces	5 pieces	5 easy to transport pieces	5 easy to transport pieces

**WARRANTY:** The Celebrity, Legend, XL, Sundancer and Sidekick Scooter Lines all have a 3 Year Limited Warranty.