

# TriShell®

1656



TriShell is comfortable, lightweight and conforms to the body and the bra for a totally natural feeling. It includes an outer silicone shell, small inner cotton filler and double thick cotton back cover. Thus it can be worn as a total breast replacement, as a partial, or as an enhancer by using only the soft silicone shell. It answers the needs of almost any breast problem, and can be worn with a feeling of total confidence day or night. Sizes 3 through 12.

# TriVaria®

1640



TriVaria's symmetrical design offers the same precise fit as a super-premium asymmetric form. The major extension has been sculpted to naturally conform to the body giving the apex area the look and feel of natural breast tissue. Accommodates a wide range of different surgeries; can be worn in any properly supporting bra. Provides more chest wall coverage, giving the wearer a perfect, firm fit. Sizes 2 through 12.

# AIRWAY®

Division of Surgical Appliance Industries Inc.  
Cincinnati, Ohio 45209  
(513) 271-4594 or 800-868-0458

In Canada, contact:  
Airway Surgical Appliances Ltd.  
Nepean, Ontario K2E 7J4  
(613) 723-4790 or 800-267-3476 (Central)  
800-267-8098 (National)

5990M02



# TriLine


An exquisite collection of specially shaped breast forms providing the ultimate in versatility and comfort

# TriTex

1655



Symmetrical 3-wing shell



- Silicone breast prosthesis (1)
- Velvet soft microfibre fabric (2)
- Transportation of moisture (3)
- Skin (4)


TriTex features velvet soft microfibre backing on the side that faces the chest wall. Together with a specially designed system of air chambers it maintains a constant flow of air, providing a pleasant, healthy contact surface with the skin.

The microfibres in TriTex guarantee previously unimagined wearing comfort. They rapidly absorb perspiration, drawing it away from the body like a wick so that it can evaporate. This lets the skin feel warm and comfortable; the clamminess of damp clothing becomes a thing of the past.

The microfibre backing is firmly bound to the prosthesis. Therefore it is recommended that TriTex be worn directly against the skin, where its comfort and quality are most evident.

The weight of TriTex has been reduced to a minimum, but its carefully sculpted contours allow for completely natural movement. It is lightweight and comfortable to wear, even in the case of lymphatic edema or heavy bosom; easy to handle and take care of.

TriTex is the gentle, skin friendly alternative that offers the client more comfort and a better quality of life!



# TriWing S

1680



A superior product at a mid-price cost. TriWing's apex along with accentuated medial and lateral flanges provide advantages of an asymmetrically shaped form. Fits the natural shape of the chest wall without depending on the pressure of a bra to create contour. Provides tissue compensating capabilities as well as a more natural nipple position. Offers the utmost in appearance, comfort and value. Sizes 0 through 10.

# TriLite

1654



TriLite is truly a new concept in a lightweight prosthesis. Its unique posterior surface design makes it especially suited for wear while swimming, as it permits unimpaired drainage upon leaving the water and thus allows for rapid drying of the skin. It provides perfectly natural contours especially in a wet or dry bathing suit, and in lightweight lingerie. Sizes 2 through 12.